

COLORIMETRIC DETERMINATION OF ASPIRIN SMART WORKSHEET

PART A. PREPARATION OF SOLUTIONS FOR BEER LAW GRAPH

COMPOSITION OF SOLUTION S

The molar mass of acetylsalicylic acid is 180.16 g mol^{-1}

Part A: Mass of acetylsalicylic acid data and molar concentration calculation

	Unrounded	Rounded	Unit
Mass of acetylsalicylic acid in 1 L solution		<div><div></div>1.6023</div>	<div><div></div>g</div>
1. Molar concentration of acetylsalicylic acid	<div><div></div>0.008893761</div>	<div><div></div>0.00889</div>	<div><div></div>mol L⁻¹</div>

PREPARATION OF DILUTED SOLUTIONS TO BE USED FOR BEER'S LAW GRAPH

Final volume is 50.00 mL for all solutions, including blank

Part A: Volume data, ferric salicylate ion concentration calculation and absorbance data

Trial	Volume of solution S (mL)	Final volume (mL)	2. [Fe-salicylate ion] (mol L ⁻¹)		Absorbance
			Unrounded	Rounded	
Blank (FeCl ₃ soln used to calibrate colorimeter)	<div><div>0.00</div></div>	<div><div>50.00</div></div>	<div><div>0.00000</div></div>	<div><div>0.00</div></div>	<div><div>0.000</div></div>
1	<div><div>2.00</div></div>	<div><div>50.00</div></div>	<div><div>0.0003557504</div></div>	<div><div>0.000356</div></div>	<div><div>0.194</div></div>
2	<div><div>3.00</div></div>	<div><div>50.00</div></div>	<div><div>0.0005336257</div></div>	<div><div>0.000534</div></div>	<div><div>0.271</div></div>
3	<div><div>4.00</div></div>	<div><div>50.00</div></div>	<div><div>0.0007115009</div></div>	<div><div>0.000712</div></div>	<div><div>0.402</div></div>
4	<div><div>5.00</div></div>	<div><div>50.00</div></div>	<div><div>0.0008893761</div></div>	<div><div>0.000889</div></div>	<div><div>0.524</div></div>
5	<div><div>6.00</div></div>	<div><div>50.00</div></div>	<div><div>0.0010672513</div></div>	<div><div>0.00107</div></div>	<div><div>0.622</div></div>

Data check

Absorbance Values: ✔ Data passes ascending order check

GRAPH CALCULATIONS

Using the graph of absorbance vs [Fe-salicylate ion], complete the table below.

Part A: Slope and intercept calculation

	Unrounded	Rounded	Unit
3. Slope	<div><div><div></div></div><div>590.542</div><div><div></div></div></div>	<div><div><div></div></div><div>5.91×10^2</div><div><div></div></div></div>	<div><div><div></div></div><div>L mol⁻¹</div><div><div></div></div></div>
4. <i>y</i> -intercept	<div><div><div></div></div><div>-0.0146429</div><div><div></div></div></div>	<div><div><div></div></div><div>-0.0146</div><div><div></div></div></div>	<div><div><div></div></div><div>Unitless</div><div><div></div></div></div>

YOUR PROGRESS ON THE PREPARATION OF SOLUTIONS FOR BEER LAW GRAPH

CORRECT	22 / 22	POINTS AWARDED98 / 98	AUTOSOLVED	0 / 22	NOT FINISHED	0 / 41
---------	---------	-----------------------	------------	--------	--------------	--------

PART B. ANALYSIS OF ASPIRIN TABLET

ASPIRIN TABLET DATA

Part B: Strength and absorbance of aspirin

	Data	Unit
Strength of aspirin tablet (mg acetylsalicylic acid)	<div><div><div></div></div><div>325</div><div><div></div></div></div>	mg
Absorbance of aspirin solution	<div><div><div></div></div><div>0.109</div><div><div></div></div></div>	Unitless

ASPIRIN TABLET CALCULATIONS

Part B: Calculations

	Unrounded	Rounded	Units
5. Molar concentration of iron salicylate ion in 50 mL volumetric flask (calculated using equation of line)	<div><div><div></div></div><div>0.000209372</div><div><div></div></div></div>	<div><div><div></div></div><div>0.000209</div><div><div></div></div></div>	mol L ⁻¹
6. Molar concentration of iron salicylate ion in 250 mL volumetric flask (parent aspirin stock solution)	<div><div><div></div></div><div>0.00523430</div><div><div></div></div></div>	<div><div><div></div></div><div>0.00523</div><div><div></div></div></div>	mol L ⁻¹
7. Moles of iron salicylate ion in 250 mL volumetric flask (parent aspirin stock solution)	<div><div><div></div></div><div>0.00130857</div><div><div></div></div></div>	<div><div><div></div></div><div>0.00131</div><div><div></div></div></div>	mol
8. Calculated mass of acetylsalicylic acid in the tablet	<div><div><div></div></div><div>235.753</div><div><div></div></div></div>	<div><div><div></div></div><div>236</div><div><div></div></div></div>	mg
9. Percent error in mass of acetylsalicylic acid in the tablet	<div><div><div></div></div><div>27.3846</div><div><div></div></div></div>	<div><div><div></div></div><div>27</div><div><div></div></div></div>	%

Quality of Data

% Error	<div><div><div></div></div><div>% error above 20%</div><div><div></div></div></div>
---------	---

YOUR PROGRESS ON

CORRECT	10 / 10	POINTS AWARDED44 / 50	AUTOSOLVED	0 / 10	NOT FINISHED	0 / 12
---------	---------	-----------------------	------------	--------	--------------	--------

YOUR OVERALL PROGRESS

Visual status toggles for statistics by question type