

ACIDS, BASES, AND TITRATION CURVES I SMART WORKSHEET

PART A: DETERMINATION OF THE CONCENTRATION OF THE NAOH TITRANT

DATA ANALYSIS

UNITS

Choose the correct units for the variables below:

Variable	Units	
Mass KHP	<div><div><div></div></div><div>g</div></div>	<div><div></div></div>
Volume NaOH used ($\Delta V = V_f - V_i$)	<div><div><div></div></div><div>mL</div></div>	<div><div></div></div>
NaOH Concentration	<div><div><div></div></div><div>mol/L</div></div>	<div><div></div></div>

DATA

Use your data collected in the lab to complete the table below.

- Trials 1 and 2 are compulsory, trials 3-5 are optional.
- The molar mass of KHP is 204.22 g/mol.
- $K_a = 1.80 \times 10^{-5}$

Trial #	Trials				
	1	2	3 (optional)	4 (optional)	5 (optional)
Mass KHP	<div><div><div></div></div><div>0.4665</div></div> <div><div></div></div>	<div><div><div></div></div><div>0.4550</div></div> <div><div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>
Volume NaOH used	<div><div><div></div></div><div>9.88</div></div> <div><div></div></div>	<div><div><div></div></div><div>7.48</div></div> <div><div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>
1. NaOH Concentration (Unrounded)	<div><div><div></div></div><div>0.231205</div></div> <div><div></div></div>	<div><div><div></div></div><div>0.297860</div></div> <div><div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>
2. NaOH Concentration (Rounded)	<div><div><div></div></div><div>0.231</div></div> <div><div></div></div>	<div><div><div></div></div><div>0.298</div></div> <div><div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>	<div><div><div></div></div></div>

MOLARITY ANALYSIS

CHOICE OF 2 TRIALS

Indicate the 2 trials used to calculate average [NaOH]:

	Choice 1	Choice 2
Trial number	<div><div><div></div></div><div>1</div></div> <div><div></div></div>	<div><div><div></div></div><div>2</div></div> <div><div></div></div>

AVERAGE MOLARITY CALCULATION

	Unrounded	Rounded
4, 5. Average Molarity of NaOH (mol/L)	<div><div><div></div></div><div>0.264500</div></div> <div><div></div></div>	<div><div><div></div></div><div>0.265</div></div> <div><div></div></div>

CORRECT

9 / 9

POINTS AWARDED

36 / 36

AUTOSOLVED

0 / 9

NOT FINISHED

0 / 15

PART B: TITRATION OF ACETIC ACID

DATA

Use your data collected in the lab to complete the table below.

Trial #	Trials		Units
	1	2	
Volume of acetic acid dispensed	8.50 ✓	8.50 ✓	mL ✓
Volume NaOH added to reach equivalence point (from 2 nd derivative analysis on LabQuest)	4.05 ✓	3.42 ✓	mL ✓
pH halfway to the equivalence point (from analysis of titration curve on LabQuest)	4.62 ✓	4.76 ✓	unitless ✓
p <i>K</i> _a of acetic acid	4.62 ✓	4.76 ✓	unitless ✓

p*K*_a ANALYSIS

AVERAGE p*K*_a

Calculate the average of measured p*K*_a values:

	Unrounded	Rounded
Average p <i>K</i> _a	4.690 ✓	4.69 ✓

CALCULATING *K*_a

Calculate the *K*_a from the average p*K*_a (4.690):

	Unrounded	Rounded
<i>K</i> _a from average p <i>K</i> _a	0.0000204174 ✓	0.000020 ✓

CALCULATING p*K*_a FROM ACCEPTED *K*_a

Calculate the p*K*_a of acetic acid from accepted *K*_a (1.80 × 10⁻⁵):

	Unrounded	Rounded
p <i>K</i> _a from accepted <i>K</i> _a	4.74473 ✓	4.74 ✓

CALCULATING % ERROR

Calculate % error of the average p*K*_a:

	Unrounded	Rounded
% error	1.15349 ✓	1.2 ✓

CONCENTRATION ANALYSIS

MOLARITY CALCULATIONS OF 2 TRIALS

Calculate the molarity of acetic acid:

Trial #	Trials		Units
	1	2	
Molarity of acetic acid (Unrounded)	0.126111 ✓	0.106463 ✓	mol/L ✓

Molarity of acetic acid (Rounded)	<div><div>0.126</div><div>✓</div></div>	<div><div>0.106</div><div>✓</div></div>	mol/L
--------------------------------------	---	---	-------

AVERAGE MOLARITY OF ACETIC ACID

Calculate the average molarity of acetic acid:

	Unrounded	Rounded
Average molarity of acetic acid	<div><div>0.116000</div><div>✓</div></div>	<div><div>0.116</div><div>✓</div></div>

% DIFFERENCE BETWEEN THE 2 TRIALS

Calculate the % difference in acetic acid between the two trials:

	Unrounded	Rounded
% difference between two trials	<div><div>17.2414</div><div>✓</div></div>	<div><div>17</div><div>✓</div></div>


YOUR PROGRESS ON PART B

<div>CORRECT</div> <div>23 / 23</div>	<div>POINTS AWARDED</div> <div>98 / 98</div>	<div>AUTOSOLVED</div> <div>0 / 23</div>	<div>NOT FINISHED</div> <div>0 / 29</div>
---------------------------------------	--	---	---

YOUR OVERALL PROGRESS

Visual status toggles for statistics by question type

☐ Units ☐ Calculation ☐ Rounding



END OF ACTIVITY