

ACIDS, BASES, AND TITRATION CURVES II SMART WORKSHEET

PART A: TITRATION OF PHOSPHORIC ACID

DATA

	Data	Unit
Average molarity of NaOH used in experiment 8	<div><div></div>0.2646<div></div></div>	<div><div></div>mol/L<div></div></div>

	Trial 1	Trial 2	Units
Volume H ₃ PO ₄ dispensed	<div><div></div>10.00<div></div></div>	<div><div></div>10.00<div></div></div>	<div><div></div>mL<div></div></div>
Volume NaOH added to reach first equivalence point (from 2 nd derivative analysis on LabQuest)	<div><div></div>4.00<div></div></div>	<div><div></div>5.30<div></div></div>	<div><div></div>mL<div></div></div>
Volume NaOH added to reach second equivalence point (from 2 nd derivative analysis on LabQuest)	<div><div></div>11.55<div></div></div>	<div><div></div>12.62<div></div></div>	<div><div></div>mL</div>
pH halfway to first equivalence point (from titration curve analysis on LabQuest)	<div><div></div>0.92<div></div></div>	<div><div></div>1.17<div></div></div>	<div><div></div>unitless<div></div></div>
pH halfway between the first and second equivalence points (from titration curve analysis on LabQuest)	<div><div></div>6.46<div></div></div>	<div><div></div>6.63<div></div></div>	<div><div></div>unitless</div>

pK_a ANALYSIS

The acid dissociation constants for phosphoric acid are:

- $K_{a1} = 0.0075$
- $K_{a2} = 6.2 \times 10^{-8}$
- $K_{a3} = 4.2 \times 10^{-13}$

	Unrounded	Rounded	Units
1. Average experimental pK _{a1} of H ₃ PO ₄	<div><div></div>1.04500<div></div></div>	<div><div></div>1.05<div></div></div>	<div><div></div>unitless<div></div></div>
2. Average experimental pK _{a2} of H ₃ PO ₄	<div><div></div>6.54500<div></div></div>	<div><div></div>6.55<div></div></div>	<div><div></div>unitless</div>
3. Actual pK _{a1} of H ₃ PO ₄	<div><div></div>2.12494<div></div></div>	<div><div></div>2.12<div></div></div>	<div><div></div>unitless</div>
4. Actual pK _{a2} of H ₃ PO ₄	<div><div></div>7.20761<div></div></div>	<div><div></div>7.21<div></div></div>	<div><div></div>unitless</div>
5. % error of your average pK _{a1}	<div><div></div>50.8221<div></div></div>	<div><div></div>50.8<div></div></div>	<div><div></div>unitless<div></div></div>
6. % error of your average pK _{a2}	<div><div></div>9.19320<div></div></div>	<div><div></div>9.2<div></div></div>	<div><div></div>unitless</div>

ANALYSIS OF CONCENTRATION OF H₃PO₄

Trial 1	Trial 2	Units
---------	---------	-------

7. Concentration of H_3PO_4 (use data for second equivalence point) <i>Unrounded</i>	<div><div><div></div></div><div>0.152807</div><div><div></div></div></div>	<div><div><div></div></div><div>0.166963</div><div><div></div></div></div>	<div><div><div></div></div><div>mol/L</div><div><div></div></div></div>
Concentration of H_3PO_4 (use data for second equivalence) <i>Rounded</i>	<div><div><div></div></div><div>0.1528</div><div><div></div></div></div>	<div><div><div></div></div><div>0.1670</div><div><div></div></div></div>	<div><div><div></div></div><div>mol/L</div><div><div></div></div></div>

	Value	Units
8. average $[\text{H}_3\text{PO}_4]$ <i>Unrounded</i>	<div><div><div></div></div><div>0.159900</div><div><div></div></div></div>	mol/L
average $[\text{H}_3\text{PO}_4]$ <i>Rounded</i>	<div><div><div></div></div><div>0.1599</div><div><div></div></div></div>	mol/L

YOUR PROGRESS ON 'PART A: TITRATION OF PHOSPHORIC ACID'

CORRECT	25 / 25	POINTS AWARDED 96 / 96	AUTOSOLVED	0 / 25	NOT FINISHED	0 / 36
----------------	----------------	-------------------------------	-------------------	---------------	---------------------	---------------

PART B: TITRATION OF AN UNKNOWN DIPROTIC ACID

Code of unknown acid:

B

	Trial 1	Trial 2	Units
Mass of unknown acid	<div><div><div></div></div><div>0.1434</div><div><div></div></div></div>	<div><div><div></div></div><div>0.1595</div><div><div></div></div></div>	<div><div><div></div></div><div>g</div><div><div></div></div></div>
Volume NaOH added to reach first equivalence point	<div><div><div></div></div><div>7.82</div><div><div></div></div></div>	<div><div><div></div></div><div>8.31</div><div><div></div></div></div>	<div><div><div></div></div><div>mL</div><div><div></div></div></div>
Volume NaOH added to reach second equivalence point	<div><div><div></div></div><div>16.51</div><div><div></div></div></div>	<div><div><div></div></div><div>17.74</div><div><div></div></div></div>	<div><div><div></div></div><div>mL</div><div><div></div></div></div>

MOLAR MASS ANALYSIS

	Trial 1	Trial 2	Units
9. Molar mass of unknown acid (use your second equivalence point data) <i>Unrounded</i>	<div><div><div></div></div><div>65.6511</div><div><div></div></div></div>	<div><div><div></div></div><div>67.9590</div><div><div></div></div></div>	<div><div><div></div></div><div>g/mol</div><div><div></div></div></div>
10. Molar mass of unknown acid (use your second equivalence point data) <i>Rounded</i>	<div><div><div></div></div><div>65.65</div><div><div></div></div></div>	<div><div><div></div></div><div>67.96</div><div><div></div></div></div>	<div><div><div></div></div><div>g/mol</div><div><div></div></div></div>

	Value	Units
11. Average molar mass of unknown acid <i>Unrounded</i>	<div><div><div></div></div><div>66.8050</div><div><div></div></div></div>	g/mol
12. Average molar mass of unknown acid <i>Rounded</i>	<div><div><div></div></div><div>66.81</div><div><div></div></div></div>	g/mol

pH DATA

	Trial 1	Trial 2	Units
pH halfway to first equivalence point (from analysis of LabQuest titration curve)	<div><div><div></div></div><div>0.89</div><div><div></div></div></div>	<div><div><div></div></div><div>0.89</div><div><div></div></div></div>	<div><div><div></div></div><div>unitless</div><div><div></div></div></div>
pH halfway between 1 st and 2 nd equivalence points (from analysis of LabQuest titration curve)	<div><div><div></div></div><div>5.54</div><div><div></div></div></div>	<div><div><div></div></div><div>5.65</div><div><div></div></div></div>	<div><div><div></div></div><div>unitless</div><div><div></div></div></div>

pK_a ANALYSIS

	Unrounded	Rounded	Units
13. Average experimental pK_{a1}	<div><div></div>0.890000<div></div></div>	<div><div></div>0.89<div></div></div>	<div><div></div>unitless<div></div></div>
14. Average experimental pK_{a2}	<div><div></div>5.59500<div></div></div>	<div><div></div>5.60<div></div></div>	<div><div></div>unitless</div>

YOUR PROGRESS ON 'PART B: TITRATION OF AN UNKNOWN DIPROTIC ACID'

CORRECT

15 / 15

POINTS AWARDED

60 / 60

AUTOSOLVED

0 / 15

NOT FINISHED

0 / 26

PART C: IDENTIFICATION OF UNKNOWN ACID

Compare experimental molar mass and pK_{a1} to information for 2 possible unknowns.

Acid	pK_{a1}	Molar mass (g mol ⁻¹)
Oxalic acid dihydrate	1.27	126.07
Maleic acid	1.92	116.1

Closest match:

maleic acid

Closest match of data	<div>Only one of your pH or molar mass values is closest to that of maleic acid</div> <div>1</div>
Correct identification	<div>Well done, you correctly identified the unknown acid</div> <div>2</div>

(Report your rounded value to 3 significant figures)

Make sure to use the literature molar mass of the correct acid in the % error calculation. The acid you selected for the 'Closest match' may not be the correct acid.

	Unrounded	Rounded
15. % error in molar mass of unknown acid	<div><div></div>42.4591<div></div></div>	<div><div></div>42.5<div></div></div>
Quality of data	<div>% difference above 10%</div> <div>0</div>	

YOUR PROGRESS ON 'PART C: IDENTIFICATION OF UNKNOWN ACID'

CORRECT

5 / 5

POINTS AWARDED

13 / 24

AUTOSOLVED

0 / 5

NOT FINISHED

0 / 3

YOUR OVERALL PROGRESS

Visual status toggles for statistics by question type

☐ Units ☐ Calculations ☐ Rounding ☐ Quality of data

